

ENVILIFE

CB-5HP-C

Digital Uniform Temperature Hot Plate



Features

- Aluminum hot plate, uniform surface temperature
- Heating up fast and uniform, heating power is stepless adjustable
- Digital hot plate has the more functions as follows
- Digital display can read the temperature more clearly

SPECIFICATIONS

• Working time	Continuous
• Max Temperature	350°C
• Structure	
• Outer Shell	Iron plate welding, surface spraying
• Hot Plate	Aluminum tensile
• Driving Mode	No
• Heating Power	1200W
• Controller	
• Continuous operation	Yes
• Setting method	Digital
• Timer	Digital
• Operation function	No
• Temperature sensor	K model
• Exterior size (W*L*H)	30×30×18 cm
• Packing Size (W*L*H)	34×37×22 cm
• Plate Size	300×300 mm
• Plate load	20Kg
• Power Supply (50/60Hz)	220V
• NW/GW	6.2/7 Kg